


GX-271 ASPEC™ with VERITY® 4060 Single Syringe Pump


Plumbing Diagram for No Air/Gas Configuration



VERITY® 4060 Single Syringe Pump Diagram

NOTE Store the excess transfer tubing in the compartment under the syringe pump.

- 1.1 mL transfer tubing (part number 499424013)
- 5 mL transfer tubing (part number 499671142)
- 10 mL transfer tubing (part number 499672282)
- 25 mL transfer tubing (part number 4996757021)



GX-271 ASPEC™ with VERITY® 4060 Single Syringe Pump

Plumbing Diagram for Air/Gas Configuration

