






Gilson Quality Control Report Pipette Head Addendum

PIPETMAX® 268 Pipette Head

Default Calibration Values

	<i>MAX8X20</i>			
	<i>Volume Setting</i>	<i>1 µL Minimum</i>	<i>10 µL 50%</i>	<i>20 µL Nominal</i>
	Calibration Volume Adjustment	1000	10000	20000
	Calibration Offset	990	8140	16140
	<i>MAX8X200</i>			
	<i>Volume Setting</i>	<i>20 µL Minimum</i>	<i>100 µL 50%</i>	<i>200 µL Nominal</i>
	Calibration Volume Adjustment	20000	100000	200000
	Calibration Offset	1660	8040	15960
	<i>MAX1000</i>			
	<i>Volume Setting</i>	<i>100 µL Minimum</i>	<i>500 µL 50%</i>	<i>1000 µL Nominal</i>
	Calibration Volume Adjustment	99500	496000	994000
	Calibration Offset	1680	8080	16064

Please refer to the *PIPETMAX® 268 User's Guide* (part number LT255519), Add Pipette Head(s) section for information about scanning or entering pipette head calibration values.

